





Nasal Implants



NASAL DORSAL Lacks a Columella

| CAT # | A | B | C |
|---------|------|-----|------|
| SND1-1F | 4.25 | 1.0 | 0.15 |
| SND1-2F | 4.3 | 1.1 | 0.25 |
| SND1-3F | 4.4 | 1.1 | 0.35 |
| SND1-4F | 4.4 | 1.1 | 0.45 |





Eliminates the columella



VOLOSHIN NASAL Dorsal Columella Design

| CAT # | A | B | C | D | E |
|-------|-----|-----|------|------|------|
| SV1 | 4.9 | 2.8 | 1.05 | 0.75 | 0.25 |
| SV2 | 5.4 | 2.8 | 1.05 | 0.75 | 0.25 |
| SV3 | 4.8 | 3.1 | 1.05 | 0.70 | 0.45 |
| SV4 | 5.5 | 3.2 | 1.1 | 0.69 | 0.50 |
| SV5 | 5.1 | 3.2 | 0.87 | 0.70 | 0.40 |






NASAL STYLES 1 THRU 16

| CAT # | A | B | C | D | E |
|--------|-----|-----|------|------|------|
| SN-1 | 5.6 | 2.3 | 0.85 | 0.65 | 0.6 |
| SN-2 | 5.2 | 2.7 | 1.0 | 0.7 | 0.55 |
| SN-3 | 4.5 | 2.7 | 0.9 | 0.6 | 0.45 |
| SN-4 | 4.9 | 1.7 | 0.93 | 0.6 | 0.35 |
| SN-5 | 4.9 | 1.9 | 1.05 | 0.75 | 0.45 |
| SN-6 | 5.2 | 1.9 | 0.85 | 0.65 | 0.45 |
| SN-7-1 | 5.5 | 2.5 | 0.95 | 0.5 | 0.35 |
| SN-7-2 | 5.5 | 2.5 | 0.95 | 0.5 | 0.25 |
| SN-8 | 4.7 | 2.5 | 1.05 | 0.65 | 0.35 |
| SN9 | 4.8 | 2.9 | 1.05 | 0.5 | 0.55 |

Hanson Medical offers the largest selection of nasal implant styles and sizes in the market. They are available in flesh tone or white. Most implants are easily modified with a scalpel. The columellar strut is almost always trimmed down, if not eliminated completely. Patient selection is critical. The patient should have adequate dermal and subcutaneous skin layers to cover the implant. Free sizers are available for each style of implant ordered.

| CAT # | A | B | C | D | E |
|--------|-----|-----|------|-----|------|
| SN10-1 | 4.9 | 3.0 | 0.95 | 0.6 | 0.45 |
| SN10-2 | 5.1 | 3.1 | 1.2 | 0.8 | 0.5 |
| SN11 | 4.6 | 2.8 | 1.1 | 0.6 | 0.5 |
| SN12-1 | 4.6 | 2.6 | 1.0 | 0.5 | 0.3 |
| SN12-2 | 5.0 | 3.1 | 1.1 | 0.7 | 0.3 |
| SN14 | 4.5 | 2.5 | 0.8 | 0.6 | 0.35 |
| SN16 | 5.0 | 2.6 | 0.9 | 0.7 | 0.3 |



Hanson Medical Inc.
PO Box 1160, Kingston, WA 98346, USA
Toll Free (800) 771-2215
Local (360) 297-1997
Fax (360) 297-1998
Email info@hansonmedical.com